

```
---
title: "Exploring FCM Data In R"
author: "R. Burke Squires, Radina Droumeva"
date: "7/11/2017"
output: ioslides_presentation
---
```

LICENSE

```
```{r setup, include=FALSE}
knitr::opts_chunk$set(echo = FALSE)
```
```

## R Markdown

This is an R Markdown presentation. Markdown is a simple formatting syntax for authoring HTML, PDF, and MS Word documents. For more details on using R Markdown see <http://rmarkdown.rstudio.com>.

When you click the **Knit** button a document will be generated that includes both content as well as the output of any embedded R code chunks within the document.

#####

# Flow Cytometry data

## Look at the help files to search for a function:

```
```{r}
?read.FCS
??read.FCS
```
```

## Load the package which extends the functionality of R to work with flow data

```
```{r}
library(flowCore)
```
```

```
## Make sure R knows which directory our data will be read from
```

```
```{r}  
getwd()  
setwd('./data')  
dir()  
dir('fullFCS/')  
```
```

```
## Read an FCS file
```

```
```{r}  
f <- read.FCS('./data/fullFCS/100715.fcs')  
f  
# 'f' is a flowFrame object. See ?flowFrame for details and to  
# see what you can do with it  
```
```

```
## Slide
```

```
```{r}  
# how many events the file has  
nrow(f)  
# the channel names:  
colnames(f)  
# Extract the expression values into a matrix  
E <- exprs(f)  
```
```

```
```{r}  
dim(E)  
```
```

```
## The expression values are like a matrix
```

```
-- Each cell has a row of measurements - one for each channel.  
Here are the first 10 cells:
```

```
```{R}  
E[1:10, ]  
```
```

```
## Explore the meta data stored within the FCS file
```

```

```{r}
f@description
names(f@description)
f@description$`TUBE NAME`
f@parameters@data
f@parameters@data[1, c("minRange", "maxRange")]
```

```

## Try a simple plot

-- note the error R gives you.  
It says that you have to first load the 'flowViz' library  
before you can plot FCM files.

```

```{r}
plot(f, c("FSC-A", "SSC-A"))
library(flowViz)
plot(f, c("FSC-A", "SSC-A"), ylim = c(0, 5000), smooth=FALSE)
# Note SSC-A is the third parameter (P3) and the meta data
tells us it is to be viewed on a LOG scale:
colnames(f)[3] # See that this is SSC-A
f@description$`P3DISPLAY`
```

```

## Read a flow set

```

```{r}
# Now read a flow set
fs <- read.flowSet(path = './data/fullFCS', pattern = ".fcs")
fs
# You can see sample names as well as the channel names
sampleNames(fs)
length(fs)
colnames(fs)
# A flowSet object is similar to a list, a list of flowFrames
fs[["100715.fcs"]]
fs[[1]]
```

```

## fsApply

```

```{r}
# Use fsApply to get cell counts for all samples
nrow(fs[[1]])
fsApply(fs, nrow)

```

```

# Use fsApply to extract the TUBE NAME keyword in all samples
fsApply(fs, function(f) f@description$`TUBE NAME`)
```

# Plotting exercise

```{r}
### Plotting exercise
#####
plot(fs[[2]], c("FSC-A", "SSC-A"), ylim = c(0, 5000), smooth =
FALSE)
# Plot the density of the forward scatter area values for the
first sample:
E <- exprs(fs[[1]])
fscValues <- E[, "FSC-A"]
fscValues[1:10]
plot(density(fscValues))
```

## plot all 3 samples on one plot

```{r}
# We can plot all 3 samples on one plot:
par (mfrow = c(3, 1)) # This creates a plot region with a
single column of 3 subplots
plot(fs[[1]], c("FSC-A", "SSC-A"), main = sampleNames(fs)[1],
ylim = c(0, 5000), smooth=FALSE)
plot(fs[[2]], c("FSC-A", "SSC-A"), main = sampleNames(fs)[1],
ylim = c(0, 5000), smooth=FALSE)
plot(fs[[3]], c("FSC-A", "SSC-A"), main = sampleNames(fs)[1],
ylim = c(0, 5000), smooth=FALSE)
```

```